


- A. Insulating cap, white ABS thermoplastic.
- B. Base plate, for skin contact, 300 series non-allergic stainless steel.
- C. 300 series stainless steel tabs secure fabric material in place.
- D. 300 series stainless steel hinge maintains electrical continuity even when open.
- E. Eccentric hinge* allows wide opening for easy adjustment and exerts significant force on band material to secure against base plate when closed.
- F. 4mm or 10mm snap stud machined from 300 series stainless steel. Fits to 4mm or 10mm spring loaded medical type socket. Coil cord disengages from snap with .45Kg - 2.27Kg of perpendicular force. Durable, providing long life well beyond that of sheet metal snaps.
- G. Knitted elastic band material of light blue fiber, 19 mm wide, 2.54 mm thick and 241 mm long. Material made up of 30 strands of spandex, 8 strands of silver coated nylon 6 yarn, and 41 strands of textured polyester.
- H. The conductive silver plated nylon 6 represents 10% of the material by weight. It is configured into 4 separate bands of 2 strands each for a total of 8 strands which are randomly knit in the inner surface (blue side).
- I. Manufacturing date code stamped on inside of plastic cap.
- J. Underwriters Laboratory listing 90P1.

Designed to satisfy wrist strap requirements of IEC 61340-5.

Minimum wrist circumference = 64mm

Maximum wrist circumference = 355mm

*U.S. Patent No. 4,662,695



Elastic Wrist Band Only

Colour	4mm Stud	10mm Stud
Blue	60030	60060



Wristband, Elastic, Adjustable

CHARLESWATER EUROPE LTD.
 UNIT 17. MILLBROOK BUSINESS PARK, SYBRON WAY
 CROWBOROUGH, EAST SUSSEX TN6 3JZ UNITED KINGDOM
 PHONE: 00 44 (0) 1892-665313, FAX: 00 44 (0) 1892-668838
 INTERNET: www.charleswater.co.uk

Drawing Number
60030.E

DATE:
06/01